# Homedics

**automatic arm  
blood pressure monitor**

**3 YEAR GUARANTEE**

TMB-1491-SN

## INTRODUCTION

### General Description

Thank you for selecting a Homedics arm blood pressure monitor. The monitor features blood pressure measurement, pulse rate measurement and the result storage. The design provides you with three years of reliable service. Readings taken by the blood pressure monitor are equivalent to those obtained by a trained observer using the cuff and stethoscope auscultation method.

This manual contains important safety and care information, and provides step by step instructions for using the product. Read the manual thoroughly before using the product.

### Features:

- 60mm x 40.5 mm LCD
- Maximum 2 users, 99 records each
- Measuring during inflation

### Indications for Use

This Blood Pressure Monitor is a digital monitor intended for use in measuring blood pressure and pulse rate with arm circumference ranging from 22cm to 32cm (about 8¾" - 12½"), 22cm to 42cm (about 8¾" - 16½"). It is intended for adult indoor use only.

### Measurement Principle

This product uses the Oscillometric Measuring method to detect blood pressure. Before every measurement, the unit establishes a "zero pressure" equivalent to the air pressure. Then it starts inflating the arm cuff, meanwhile, the unit detects pressure oscillations generated by beat-to-beat pulsatile, which is used to determine the systolic and diastolic pressure, and also pulse rate.

### Safety Information

The below signs might be in the user manual, labeling or other components. They are the requirement of standard & using.

|  |                                                                                                                                                                                                                                                                           |  |                                                                                                                                                                                                                                                            |
|--|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Refer to instruction manual/booklet<br>To signify that the instruction manual/<br>booklet must be read.                                                                                                                                                                   |  | Type BF applied part                                                                                                                                                                                                                                       |
|  | CE marking indicates that a product has been<br>assessed by the manufacturer and deemed to<br>meet EU safety, health and environmental<br>protection requirements. It is required for<br>products manufactured anywhere in the<br>world that are then marketed in the EU. |  | The symbol indicates that the product should<br>not be discarded as unsorted waste but must<br>be sent to separate collection facilities for<br>recovery and recycling.                                                                                    |
|  | Symbol for "MANUFACTURER"                                                                                                                                                                                                                                                 |  |                                                                                                                                                                                                                                                            |
|  | Symbol for "SERIAL NUMBER"                                                                                                                                                                                                                                                |  | Authorized Representative in the European<br>Community                                                                                                                                                                                                     |
|  | Symbol for "DIRECT CURRENT"                                                                                                                                                                                                                                               |  | Symbol for "Recycle"                                                                                                                                                                                                                                       |
|  | Symbol for "MANUFACTURE DATE"                                                                                                                                                                                                                                             |  | Caution: Indicates that caution is necessary<br>when operating the device or control close<br>to where the symbol is placed, or that the<br>current situation needs operator awareness<br>or operator action in order to avoid<br>undesirable consequences |
|  | Batch code                                                                                                                                                                                                                                                                |  | Symbol for "MEDICAL DEVICE"                                                                                                                                                                                                                                |
|  |                                                                                                                                                                                                                                                                           |  | Recycling advice for FR, IT and UK<br>Consumers                                                                                                                                                                                                            |

## PRECAUTION

- This device is intended for indoor, home use and is not intended for self-use in public areas.
- This device is portable, but it is not intended for use during patient transport.
- This device is not suitable for continuous monitoring during medical emergencies or operations.
- This device is intended for non-invasive measuring and monitoring of arterial blood pressure.

It is not intended for use on extremities other than the arm, or for any purpose other than obtaining a blood pressure measurement.

- This device is for adults. Do not use this device on neonates or infants. Do not use it on children and adolescents unless otherwise instructed by a medical professional.
- Consult with your physician before using this monitor if you suffer from the following conditions: common arrhythmias such as premature ventricular beats or atrial fibrillation; peripheral arterial disease; pregnancy; preeclampsia; implantation with electrical devices; undergoing intravascular therapy; arteriovenous shunt or mastectomy.

Please note that any of these conditions may affect measurement readings, in addition to patient motion, trembling or shivering.

- Do not use this device for diagnosis or treatment of any health problems or diseases. Please consult with your physician first whether the blood pressure or pulse rate readings can be used as an input in determining clinical actions. Please note that clinical actions can only be determined by the physician, otherwise it may lead to delayed treatment or other dangerous situations.
- If you are taking medication, consult your physician to determine the proper time to measure your blood pressure.
- This device may be used only for the intended use described in this manual, the manufacturer shall have no liability for any incidental, consequential, or special damages caused by misuse or abuse.
- Please use the device under the environment which is provided in the user manual. Otherwise, the performance and lifetime of the device will be impacted and reduced.
- The device may require up to 30 minutes to warm up / cool down from the minimum/ maximum storage temperature before it is ready for use.
- The blood pressure monitor and the cuff are suitable for use within the patient environment.
- Do not wash the cuff in a washing machine or dishwasher!

- The device contains sensitive electronic components. To avoid measurement errors, avoid taking blood pressure measurements near a strong electromagnetic field radiated interference signal or electrical fast transient/burst signal.
- Blood Pressure Monitor is intended for use by medical staffs and lay persons, and patient is also an intended user or operator.

## CAUTION

- Do not attempt to repair the unit yourself if it malfunctions. Only have repairs carried out by authorized service centers.
- It is recommended that the performance should be checked after repair, maintenance, and every two years of use, by retesting the requirements in limits of the error of the cuff pressure indication and air leakage (testing at least at 50 mmHg and 200 mmHg). Please contact manufacturer or distributor for authorized service personnel.
- Store your device, cuff in a clean and dry place, protect it against extreme moisture, heat, lint, dust and direct sunlight. Never place any heavy objects on it.
- Make sure the rubber tube of the cuff is not squeezed, stretched, or kinked during storage.
- Dispose of accessories, detachable parts, and the device according to the local guidelines.

## WARNING

Do not apply the cuff on an arm that has an intravenous drip or a blood transfusion attached.

- Do not kink, fold, stretch, compress, or otherwise deform the tube during measuring, as the cuff pressure might continuously increase, which could prevent blood flow and result injury.
- Taking blood pressure measurements too frequently could disrupt blood circulation and cause injuries.
- Do not apply cuff to areas on patient where skin is delicate or damaged. Check cuff site frequently for irritation.
- Do not place the cuff on the arm of a person whose arteries or veins are undergoing medical treatment, i.e. intra-vascular access or intra-vascular therapy or an arteriovenous (A-V) shunt, which could disrupt blood circulation and cause injuries.
- Do not place the cuff on the arm on the same side of a mastectomy (especially when lymph nodes have been removed). It is recommended to take measurements on the unaffected side.
- Do not wrap the cuff on the same arm to which another

monitoring device is applied. One or both devices could temporarily stop functioning if you try to use them at the same time.

- Please check (for example, by observation of the limb concerned) that the operation of the device does not result in prolonged impairment of patient blood circulation.
- On the rare occasion of a fault causing the cuff to remain fully inflated during measurement, loosen and remove the cuff immediately. Prolonged high pressure applied to the arm (cuff pressure >300 mmHg or constant pressure >15 mmHg for more than 3 minutes) might lead to bruising and discolored skin.
- Do not use this device with high-frequency (HF) surgical equipment at the same time.
- This device is not used in conjunction with oxygen rich environments, not intended for use with flammable anaesthetics, not intended for use in conjunction with flammable agents.
- Do not touch output of the batteries and the user simultaneously.
- The power cord is considered the disconnect device for isolating this equipment from supply mains.

Do not position the equipment so that it is difficult to reach or disconnect.


- Do not use this device if you are allergic to polyester, nylon, or plastic.
- Only use accessories approved by manufacturer. Using unapproved accessories might cause damage to the unit and injure users.
- If you experience discomfort during a measurement, such as pain in the arm or other complaints, press the Power button immediately to release the air from the cuff.
- Do not use the device while under maintenance, or being serviced.
- The air tube poses a risk of strangulation. Furthermore, the small parts of product and batteries present a choking hazard if swallowed. They should therefore always be kept away from infants/children.
- Sensor degradation or looseness may reduce performance of device or cause other problems.

frequency of washing, skin condition, and storage state.

- Please report to the manufacturer and the competent authority of the Member State / the FDA in which you are established about any serious incident that has occurred in relation to this device.

## NOTICE

- At the request of authorized service personnel, circuit diagrams, component part lists, descriptions, and calibration procedures will be made available by the manufacturer or distributor.
- The expected lifetime of the cuff may vary by the


| SYMBOL                                | DESCRIPTION                    | EXPLANATION                                                                                                                                                                                                             |
|---------------------------------------|--------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>SYS</b>                            | Systolic blood pressure        | High blood pressure                                                                                                                                                                                                     |
| <b>DIA</b>                            | Diastolic blood pressure       | Low blood pressure                                                                                                                                                                                                      |
| <b>Pul/min</b>                        | Pulse display                  | Pulse in beats per minute                                                                                                                                                                                               |
| <b>▼</b>                              | Deflation symbol               | The cuff is deflating.                                                                                                                                                                                                  |
| <b>M 88</b>                           | Memory                         | Indicate it is in the memory mode and which group of memory it is.                                                                                                                                                      |
| <b>AVG</b>                            | Average value                  | The average of last 3 measurements taken within 30 minutes                                                                                                                                                              |
| <b>mmHg</b>                           | mmHg                           | Measurement Unit of the blood pressure                                                                                                                                                                                  |
| <b>BAT Lo+</b>                        | Low battery                    | Batteries are low and need to be replaced                                                                                                                                                                               |
| <b>Irregular pulse rate</b>           | Irregular pulse rate           | Blood pressure monitor is detecting an irregular pulse rate during measurement.                                                                                                                                         |
| <b>Pulse rate</b>                     | Pulse rate                     | Blood pressure monitor is detecting a pulse rate during measurement.                                                                                                                                                    |
| <b>00:00 AM</b>                       | Current Time                   | Year/Month/Day, Hour/Minute                                                                                                                                                                                             |
| <b>1 2</b>                            | User 1/User 2                  | Start measurement for User 1/User 2                                                                                                                                                                                     |
| <b>Excessive Body Motion Detector</b> | Excessive Body Motion Detector | Appears when talking, moving, or shaking of the arm with the cuff on is detected during a measurement.<br>NOTE: The measured blood pressure reading may not be accurate when this symbol is displayed with the reading. |

## Monitor Components

### Component list of pressure measuring system

- 1 Cuff
- 2 Air pipe
- 3 PCBA
- 4 Pump
- 5 Valve


### List

1. Automatic Arm Blood pressure monitor


3. 4×AAA batteries

Replace the batteries under the following circumstances:


- +Lo displays on the LCD
- The LCD display is dim
- When powering on the monitor, the LCD doesn't light up.

2. Cuff (Type BF applied part) 22cm-32cm)  
(22-42cm XL cuff is available to purchase on the Homedics website)

(Please use HoMedics Authorized cuff. The size of the actual cuff please refer to the label on the attached cuff.)

4. User manual

### Installing and replacing the batteries

- Open the battery cover.
- Install the batteries as indicated in the battery compartment.  
(Always select the authorized / specified battery: 4 x AAA-size batteries).
- Replace the battery cover.


### CAUTION

- Do not use new and used batteries together.
- Do not use different types of batteries together.
- Do not dispose the batteries in fire. Batteries may explode or leak.
- Remove batteries if the device is not likely to be used for some time.
- Worn batteries are harmful to the environment. Do not dispose with daily garbage.
- Remove the old batteries from the device following your local recycling guidelines.


## SETTING DATE, TIME AND MEASUREMENT UNIT

It is important to set the clock before using your blood pressure monitor, so that a time stamp can be assigned to each record that is stored in the memory. (The setting range of the year :2021–2050 time format: 12H/24H)


1. When the monitor is off, long press "SET" button, the display will show date format "D M" or "M D". (M for month, D for day.) Press the "MEM" or "SET" button to toggle the date format between "D M" and "M D". Confirm the entry by pressing the "START/STOP" button.

2. When the display show a blinking number representing the [YEAR], change the [YEAR] by pressing the "MEM" or "SET" button. Each press of the "MEM" or "SET" button will increase or decrease the number by one in a cycling manner.


3. When you get the right year, press "START/STOP" button to confirm the entry. The screen will then show a blinking number representing the [MONTH].


4. Repeat steps 2 and 3 to set the [MONTH] and [DAY].

5. Repeat steps 2 and 3 to confirm the timeformat between [12H] and [24H].

6. Repeat steps 2 and 3 to set the [HOUR] and [MINUTE].


7. After the [MINUTE] is set, the LCD will display "done", then it will power off after several seconds.

## BEFORE YOU START

### Select the User

- When the monitor is off, press and hold "MEM" button, the user ID will show. Then press "MEM" or "SET" button to toggle the user ID between user 1 and user 2.



- Press "STRAT/STOP" button to ensure your choice, it will show User ID and "done" and then power off.

## MEASUREMENT

### Tie the Cuff

- Remove all jewelry, such as watches and bracelets from your left arm. Note: If your doctor has diagnosed you with poor circulation in your left arm, use your right arm.
- Roll or push up your sleeve to expose the skin. Make sure your sleeve is not too tight.
- Hold your arm with your palm facing up and tie the cuff on your upper arm, then position the tube off-center toward the inner side of arm in line with the little finger. Or position the artery mark  $\phi$  over the main artery (on the inside of your arm). Note: Locate the main artery by pressing with 2 fingers approximately 2 cm above the bend of your elbow on the inside of your left arm. Identify where the pulse can be felt the strongest. This is your main artery.
- The cuff should be snug but not too tight. You should be able to insert one finger between the cuff and your arm.

- Sit comfortably with your tested arm resting on a flat surface. Place your elbow on a table so that the cuff is at the same level as your heart. Turn your palm upwards. Sit upright in a chair, and take 5-6 deep breaths.
- Helpful tips for Patients, especially for patients with Hypertension:

- Rest for 5 minutes before first measuring.
- Wait at least 3 minutes between measurements. This allows your blood circulation to recover.
- Take the measurement in a silent room.
- The patient must relax as much as possible and do not move and talk during the measurement procedure.
- The cuff should maintain at the same level as the right atrium of the heart.
- Please sit comfortably. Do not cross your legs and keep your feet flat on the ground.
- Keep your back against the backrest of the chair.
- For a meaningful comparison, try to measure under similar conditions.

For example, take daily measurements at approximately the same time, on the same arm, or as directed by a physician.


### Start the Measurement


- When the monitor is off, press the "START/STOP" to turn on the monitor, and it will finish the whole measurement. (Take User 1 for example)

**START  
STOP**


LCD display


Adjust the zero

Inflating and  
measuring

Display and save  
the results.


2. Press the "START/STOP" to power off, otherwise it will turn off within 1 minute.

START  
STOP

**Note:** during the measurement, If the irregular pulse rate is detected, the symbol  will display in the measurement result.

During the measurement, once the excessive body motion, especially of the arm the cuff is worn is detected, the symbol  will flash about 5 seconds and detect again. If not detected no longer, the symbol will disappear; If still detected, the Excessive Body Motion Detector symbol will final display in the measurement result.


Tips: The measured blood pressure reading may not be accurate if this symbol is displayed.

Tips: Any time if you want to stop the measurement, press "START/STOP" button.


## DATA MANAGEMENT

### Recalling Records

1. When the monitor is off, press "MEM" button, the current user ID will display and blink.  
At this time you can toggle the user between the User 1 and User 2 by pressing "MEM" or "SET" button.  
Press "START/STOP" button to confirm and see the average of last 3 measurements.

2. Press the "MEM" or "SET" to get the record you want.

The date and time of the record will be shown alternately

The current no. is no.2  
The corresponding date is January 1st  
The corresponding time is pm 10:08

### CAUTION

Both User 1 and User 2 can store a maximum of 99 records each

## Deleting records


Delete the records by following the steps below.

(Example shown below for User 1)

1. When the monitor is in the memory recall mode, press and hold "MEM" button for about 3 seconds, the LCD will display a blinking "dEL ALL".

2. Press "START/STOP" button to confirm the deletion and the LCD will display "do nE". Once deleted, your readings cannot be restored. The LCD will display "----" like the image below.

Before you confirm the deletion, you could press "MEM" or "SET" button to switch the LCD display from "dEL ALL" to "dEL no" and then press "START/STOP" button to stop clearing the memory.

4. Press "START/STOP" button to turn off the monitor, otherwise it will power off automatically after about 1 minute.

## INFORMATION FOR USER

### Tips for Measurement

Measurements may be inaccurate if taken in the following circumstances.

Within 1 hour after eating or drinking

Immediate measurement after tea, coffee, smoking

Within 20 minutes after taking a bath

When talking or moving your fingers

In a very cold environment

When you want to discharge urine

### Maintenance

In order to get the best performance, please follow the instructions below.

#### 1. Cleaning Process:

Step 1: Make sure to switch off and unplug the device prior to cleaning.

Step 2: Use a soft, damp cloth with soapy water to clean the cuff first, and then use a soft cloth damp with clear water to remove residual soap until there is no visible residual contaminants. Attention should be paid to avoid liquid soaking into the cuff.

Step 3: Use a dry soft cloth to wipe the cuff, in order to remove residual moisture.

Step 4: Dry the cuff in a well-ventilated place after cleaning.

## 2. Disinfection Process:

- Step 1: Make sure to switch off and unplug the device prior to disinfection.
- Step 2: Use a soft, damp cloth with 70% isopropanol to disinfect the cuff for about 3 minutes. Attention should be paid to avoid liquid soaking into the cuff.
- Step 3: Use a clean dry cloth or towel to wipe off the disinfectant until there is no visible residue.
- Step 4: Dry the cuff in a well-ventilated place after disinfection.

## Suggestion:

Frequency of Cleaning and Disinfection:

For single patient multiple use, it's recommended to clean the device surface once a month or whenever it's necessary.

For multiple patient multiple use, it's recommended to clean the device every time before and after usage. Maintenance procedures shall be taken as per instruction.

## ABOUT BLOOD PRESSURE

### What are systolic pressure and diastolic pressure?

When ventricles contract and pump blood out of the heart, the blood pressure reaches its maximum value in the cycle, which is called systolic pressure. When the ventricles relax, the blood pressure reaches its minimum value in the cycle, which is called diastolic pressure.



### What is the standard blood pressure classification?

The blood pressure classification published by World Health Organization (WHO) and International Society of Hypertension (ISH) in 1999 is as follows:

#### CAUTION

Only a physician can tell your normal BP range. Please contact a physician if your measuring result falls out of the range. Kindly note that only a physician could tell whether your blood pressure value has reached a dangerous point.

| Blood pressure (mm Hg) \ Level | Optimal | Normal  | High-normal | Mild    | Moderate | Severe |
|--------------------------------|---------|---------|-------------|---------|----------|--------|
| SYS                            | <120    | 120-129 | 130-139     | 140-159 | 160-179  | ≥180   |
| DIA                            | <80     | 80-84   | 85-89       | 90-99   | 100-109  | ≥110   |

### Irregular Pulse Rate Detector

An irregular pulse rate will be detected if there is an irregular pulse rhythm while measuring systolic and diastolic blood pressure. When measurements were performed, the monitor will record all pulse intervals and calculate the average. If two or more pulse intervals were recorded, and the difference between each interval and the average is larger than  $\pm 25\%$  of the average; or if four or more pulse intervals were recorded, and the difference between each interval and the average is larger than  $\pm 15\%$  of the average value, the irregular pulse symbol will be displayed along with measurement results.

## CAUTION

The appearance of the IPR icon indicates that a pulse irregularity consistent with an irregular pulse rate was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the irregular pulse rate detector results cannot be used directly for clinical judgement.

Please seek medical advice from professionals before making any medical decisions.

### Why does my blood pressure fluctuate throughout the day?

1. Individual blood pressure varies multiple times everyday. It is also affected by the way you tie your cuff and your measurement position, so please take the measurement under the same conditions.
2. If the person takes medication, the pressure will vary more.
3. Wait at least 3 minutes for another measurement.

### Why do I get a different blood pressure at home compared to the hospital?

The blood pressure is different even throughout the day due to weather, emotion, exercise etc. Also, there is the "white coat" effect, which means blood pressure usually increases in clinical settings.

### Is the result the same if measuring on the right arm?

It is ok for both arms, but there will be some differences for some people.


We suggest you measure the same arm every time.

What you need to pay attention to when you measure your blood pressure at home:

- If the cuff is tied properly.
- If the cuff is too tight or too loose.
- If the cuff is tied on the upper arm.
- If you feel anxious.

Taking 2-3 deep breaths before beginning will be better for measuring.

Advice: Relax yourself for 4-5 minutes until you calm down.

## TROUBLESHOOTING

This section includes a list of error messages and frequently asked questions for problems you may encounter with your blood pressure monitor. If the products is not operating as you think it should, check here before arranging for servicing.

| PROBLEM         | SYMPTOM                                                                                                      | CHECK THIS                                                                     | REMEDY                                                                                                                                                                                                  |
|-----------------|--------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| No power        | Display will not light up.                                                                                   | Batteries are exhausted.                                                       | Replace with new batteries                                                                                                                                                                              |
|                 |                                                                                                              | Batteries are inserted incorrectly.                                            | Insert the batteries correctly                                                                                                                                                                          |
| Low batteries   | Display is dim or show  +Lo | Batteries are low.                                                             | Replace with new batteries                                                                                                                                                                              |
| Error message   | E 1 shows                                                                                                    | The cuff is too tight or too loose.                                            | Readjust the cuff ,not too loose or too tight and then measure again.                                                                                                                                   |
|                 | E 2 shows                                                                                                    | The monitor detected motion, talking or the pluse is too poor while measuring. | Relax for a moment and then measure again.                                                                                                                                                              |
|                 | E 3 shows                                                                                                    | The measurement process does not detect the pulse signal.                      | Loosen the clothing on the arm and then measure again                                                                                                                                                   |
|                 | E 4 shows                                                                                                    | The treatment of the measurement failed.                                       | Relax for a moment and then measure again.                                                                                                                                                              |
|                 | EExx, shows on the display.                                                                                  | A calibration error occurred.                                                  | Retake the measurement. If the problem persists, contact the retailer or our customer service department for further assistance. Refer to the warranty for contact information and return instructions. |
| Warning message | "out " shows                                                                                                 | Out of measurement range                                                       | Relax for a moment message measure again. If the problem persists, contact your physician.                                                                                                              |

## SPECIFICATIONS

|                                        |                                                                                                                                                                                                                                   |
|----------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Power supply                           | 6VDC 4×AAA batteries                                                                                                                                                                                                              |
| Display mode                           | LCD V.A.60mm×40.5mm                                                                                                                                                                                                               |
| Measurement mode                       | Oscillographic testing mode                                                                                                                                                                                                       |
| Measurement range                      | Rated cuff pressure:<br>0mmHg~299mmHg (0kPa ~ 39.9kPa)<br>Measurement pressure:<br>SYS: 60mmHg~230mmHg (8.0kPa~30.7kPa)<br>DIA: 40mmHg~130mmHg (5.3kPa~17.3kPa)<br>Pulse value: (40~199)beat/minute                               |
| Accuracy                               | Static Pressure: 5°C~40°C within ±3mmHg<br>Pulse value: ±5%<br>Clinical validation: Mean difference within ±5mmHg;<br>Standard deviation ≤8mmHg                                                                                   |
| Normal working condition               | A temperature range of :+5°C to +40°C<br>A relative humidity range of 15% to 90%, non-condensing, but not requiring a water vapour partial pressure greater than 50 hPa<br>An atmospheric pressure range of : 700 hPa to 1060 hPa |
| Storage & transportation condition     | Temperature:-20°C to +60°C<br>A relative humidity range of ≤ 93%, non-condensing, at a water vapour pressure up to 50hPa<br>An atmospheric pressure range of 500 hPa to 1060 hPa                                                  |
| Measurement perimeter of the upper arm | About 22cm~32cm                                                                                                                                                                                                                   |
| Net Weight                             | Approx.169g (Excluding the dry cells)                                                                                                                                                                                             |
| External dimensions                    | Approx.110mm×110mm×41mm                                                                                                                                                                                                           |
| Attachment                             | 4×AAA batteries,user manual                                                                                                                                                                                                       |
| Mode of operation                      | Continuous operation                                                                                                                                                                                                              |
| Degree of protection                   | Type BF applied part                                                                                                                                                                                                              |
| Protection against ingress of water    | IP21 It means the device could be protected against solid foreign objects of 12,5mm Φ and greater, and against vertically falling water drops.                                                                                    |
| Device Classification                  | Battery Powered Mode: Internally Powered ME Equipment                                                                                                                                                                             |
| Software version                       | A01                                                                                                                                                                                                                               |
| Expected Lifetime                      | Device: 3 years or 10,000 measurements (may vary based on usage conditions)<br>Cuff: 10000 times<br>Alkaline battery: About 200~300 times                                                                                         |
| Types of use/reuse                     | Multiple patient multiple use                                                                                                                                                                                                     |

WARNING: No modification of this equipment is allowed.

## CONTACT INFORMATION

For more information about our products, please visit [www.homedics.co.uk](http://www.homedics.co.uk)

Model: TMB-1491-SN

Manufactured by: Guangdong Transtek Medical Electronics Co., Ltd.

Zone A, No.105, Dongli Road, Torch Development District, 528437 Zhongshan, Guangdong, China

**Authorized European Representative:**

Medical Device Safety Service GmbH

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30175 Hannover,

Germany

**Imported into the UK by:**

FKA Brands Ltd, Somerhill Business Park, Tonbridge, Kent, TN11 0GP

**Imported into the EU by:**

FKA Brands Ltd, 29 Earlsfort Terrace, Dublin D02AY28, IE

Customer Support: [support@homedics.co.uk](mailto:support@homedics.co.uk)

## EMC GUIDANCE

The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments.

Essential performance:

Accuracy of measuring blood pressure and pulse rate

|                     |                                                                                                 |
|---------------------|-------------------------------------------------------------------------------------------------|
| Measurement Range   | Systolic pressure: 60-230 mmHg<br>Diastolic pressure: 40-130 mmHg<br>Pulse: 40-199 beats/minute |
| Rated Cuff Pressure | 0-299 mmHg (0-39.9 kPa)                                                                         |
| Accuracy            | Pressure: $\pm 3$ mmHg<br>Pulse : $\pm 5\%$                                                     |

The Basis Safety of the Blood Pressure Monitor (TMB-1491-SN) is as following:

Deviation from normal operation that poses an unacceptable risk to the patient or operator.

Warning: Don't be near the active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Table 2

| Guidance and manufacturer's declaration – electromagnetic emission                                     |                                                   |                                                   |
|--------------------------------------------------------------------------------------------------------|---------------------------------------------------|---------------------------------------------------|
| Immunity test                                                                                          | IEC 60601-1-2<br>Test level                       | Compliance level                                  |
| Electrostatic discharge (ESD)<br>IEC 61000-4-2                                                         | ±8 kV contact<br>±2 kV, ±4 kV, ±8 kV, ±15 kV, air | ±8 kV contact<br>±2 kV, ±4 kV, ±8 kV, ±15 kV, air |
| Electrical fast transient/burst<br>IEC 61000-4-4                                                       | Not applicable                                    | Not applicable                                    |
| Surge<br>IEC 61000-4-5                                                                                 | Not applicable                                    | Not applicable                                    |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>IEC 61000-4-11 | Not applicable                                    | Not applicable                                    |
| Power frequency magnetic field<br>IEC 61000-4-8                                                        | 30 A/m<br>50Hz/60Hz                               | 30 A/m<br>50Hz/60Hz                               |
| Conducted RF<br>IEC61000-4-6                                                                           | Not applicable                                    | Not applicable                                    |
| Radiated RF<br>IEC61000-4-3                                                                            | 10 V/m<br>80 MHz – 2,7 GHz<br>80 % AM at 1 kHz    | 10 V/m<br>80 MHz – 2,7 GHz<br>80 % AM at 1 kHz    |
| NOTE U <sub>T</sub> is the a.c. mains voltage prior to application of the test level.                  |                                                   |                                                   |

Table 3

| Guidance and manufacturer's declaration – electromagnetic immunity                                                       |                      |               |                                                                |                                      |                   |              |                                |                        |
|--------------------------------------------------------------------------------------------------------------------------|----------------------|---------------|----------------------------------------------------------------|--------------------------------------|-------------------|--------------|--------------------------------|------------------------|
| Radiated RF<br>IEC61000-4-3<br>(Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment) | Test frequency (MHz) | Band a) (MHz) | Service a)                                                     | Modulation b)                        | Maximum Power (W) | Distance (m) | IEC 60601-1-2 Test Level (V/m) | Compliance level (V/m) |
|                                                                                                                          | 385                  | 380-390       | TETRA 400                                                      | Pulse modulation b) 18Hz             | 1.8               | 0.3          | 27                             | 27                     |
|                                                                                                                          | 450                  | 430-470       | GMRS 460, FRS 460                                              | FM c) $\pm$ 5kHz deviation 1kHz sine | 2                 | 0.3          | 28                             | 28                     |
|                                                                                                                          | 710                  | 704-787       | LTE Band 13,17                                                 | Pulse modulation b) 217Hz            | 0.2               | 0.3          | 9                              | 9                      |
|                                                                                                                          | 745                  |               |                                                                |                                      |                   |              |                                |                        |
|                                                                                                                          | 780                  |               |                                                                |                                      |                   |              |                                |                        |
|                                                                                                                          | 810                  | 800-960       | GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5         | Pulse modulation b) 18 Hz            | 2                 | 0.3          | 28                             | 28                     |
|                                                                                                                          | 870                  |               |                                                                |                                      |                   |              |                                |                        |
|                                                                                                                          | 930                  |               |                                                                |                                      |                   |              |                                |                        |
|                                                                                                                          | 1720                 | 1700-1990     | GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4,25; UMTS | Pulse modulation b) 217 Hz           | 2                 | 0.3          | 28                             | 28                     |
|                                                                                                                          | 1845                 |               |                                                                |                                      |                   |              |                                |                        |
|                                                                                                                          | 1970                 |               |                                                                |                                      |                   |              |                                |                        |
|                                                                                                                          | 2450                 | 2400-2570     | Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7           | Pulse modulation 217 Hz              | 2                 | 0.3          | 28                             | 28                     |
|                                                                                                                          | 5240                 | 5100-5800     | WLAN 802.11 a/n                                                | Pulse modulation 217 Hz              | 0.2               | 0.3          | 9                              | 9                      |
|                                                                                                                          | 5500                 |               |                                                                |                                      |                   |              |                                |                        |
|                                                                                                                          | 5785                 |               |                                                                |                                      |                   |              |                                |                        |

### 3 YEAR GUARANTEE

FKA Brands Ltd guarantees this product from defect in material and workmanship for a period of 3 years from the date of purchase, except as noted below. This FKA Brands Ltd product guarantee does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorised accessory; alteration to the product; or any other conditions whatsoever that are beyond the control of FKA Brands Ltd. This guarantee is effective only if the product is purchased and operated in the UK / EU. A product that requires modification or adaptation to enable it to operate in any country other than the country for which it was designed, manufactured, approved and / or authorised, or repair of products damaged by these modifications is not covered under this guarantee. FKA Brands Ltd shall not be responsible for any type of incidental, consequential or special damages.

To obtain guarantee service on your product, return the product post-paid to your local service centre along with your dated sales receipt (as proof of purchase). Upon receipt, FKA Brands Ltd will repair or replace, as appropriate, your product and return it to you, post-paid. Guarantee is solely through HoMedics Service Centre. Service of this product by anyone other than the HoMedics Service Centre voids the guarantee. This guarantee does not affect your statutory rights.

For your local HoMedics Service Centre, go to [www.homedics.co.uk/servicecentres](http://www.homedics.co.uk/servicecentres)

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