HOMEDICS®

## FINGERTIP PULSE OXIMETER


## General description

Oxygen binds to hemoglobin in red blood cells when moving through the lungs. It is transported throughout the body as arterial blood. A pulse oximeter uses two frequencies of light (red and infrared) to determine the percentage (%) of hemoglobin in the blood that is saturated with oxygen. The percentage is called blood oxygen saturation, or SpO<sub>2</sub>. A pulse oximeter also measures and displays the pulse rate at the same time it measures the SpO<sub>2</sub>.

## Principle of operation principle

The principle of operation is as follows: The pulse oximeter works by applying a sensor to a pulsating arterial vascular bed. The sensor contains a dual light source and photo detector. The one wavelength of light source is 660nm, which is red light; the other is 905nm, which is infrared-red light. Skin, bone, tissue and venous vessels normally absorb a constant amount of light over time. The photo detector in finger sensor collects and converts the light into electronic signal which is proportional to the light intensity. The arterial bed normally pulsates and absorbs variable amounts of light during systole and diastole, as blood volume increases and decreases. The ratio of light absorbed at systole and diastole is translated into an oxygen saturation measurement. This measurement is referred to as SpO<sub>2</sub>.

## Diagram of operation principle

## Precautions for use

- 1. Before use, carefully read the manual.
- 2. Operation of the fingertip pulse oximeter may be affected by the use of an electro-surgical unit (ESU).
- 3. The fingertip pulse oximeter must be able to measure the pulse properly to obtain an accurate SpO<sub>2</sub> measurement. Verify that nothing is hindering the pulse measurement before relying on the SpO<sub>2</sub> measurement.
- 4. Do not use the fingertip pulse oximeter in an MRI or CT environment.
- 5. Do not use the fingertip pulse oximeter where there are restrictions on movement.
- 6. Do not use the fingertip pulse oximeter in an explosive atmosphere.
- 7. The fingertip pulse oximeter is intended only as an adjunct in patient assessment. It must be used in conjunction with other methods of assessing clinical signs and symptoms.
- 8. In order to ensure correct sensor alignment and skin integrity, the maximum application time at a single site for our device should be less than half an hour.
- 9. Do not sterilize the device using autoclaving, ethylene oxide sterilizing, or immersing the device in liquid. The device is not intended for sterilization.
- 10. Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.
- 11. This equipment complies with IEC 60601-1-2-2014 for electromagnetic compatibility for medical electrical equipment and/or systems. However, because of the proliferation of radio-frequency transmitting equipment and other sources of electrical noise in healthcare and other environments, it is possible that high levels of such interference due to close proximity or strength of a source might disrupt the performance of this device.
- 12. Portable and mobile RF communications equipment can affect medical electrical equipment
- 13. This equipment is not intended for use during patient transport outside the healthcare facility.
- 14. This equipment is not to be used adjacent to or stacked with other equipment.
- 15. It may be unsafe to:
  - use accessories, detachable parts and materials not described in the instructions for use
  - interconnect this equipment with other equipment not described in the instructions for use
  - disassemble, repair or modify the equipment
- 16. These materials that contact with the patient's skin contain medical silicone and ABS plastic enclosure are all pass the ISO10993-5 Tests for invitro cytotoxicity and ISO10993-10 Tests for irritation and delayed-type hypersensitivity.
- 17. When the signal is not readable, please do not measure.

## Contraindication

- 1. Not suitable for continuous monitoring.
- 2. According to the type of protection against electric shock: INTERNALLY POWERED EQUIPMENT;
- 3. According to the degree of protection against electric shock: TYPE BF APPLIED PART (The application part is rubber inside of the Pulse Oximeter);
- 4. According to the degree of protection against ingress of dust and water: IP22
- 5. According to the mode of operation: CONTINUOUS OPERATION

## Clinical study summary

The following details are provided to disclose actual performance observed in the clinical validation study of healthy adult volunteers. The ARMS value analysis statement and Bland-Altman plot of data is shown as following:

## Specifications

- Display Type: OLED display
- SpO<sub>2</sub>
  - Display range: 0%-100%
  - Measurement range: 70%-100%
  - Accuracy: 70%-100%±2%, 0%-69% no definition
  - Resolution: 1%
- Pulse
  - Display range: 30bpm-250bpm
  - Measure range: 30bpm-250bpm
  - Accuracy: 30bpm-99bpm, ±2bpm; 100bpm-250bpm, ±2%
  - Resolution: 1bpm
- Perfusion index
  - Display range: 0.1%-20%
  - Measure range: 0.2%-20%
  - Accuracy: 0.2%-1.0%, ±2.2digits; 1.1%-20.0%, ±20%
  - Resolution: 0.1%
- Probe LED specifications
 

	Wavelength	Radiant Power
RED	660±3nm	3.2mw
IR	905±10nm	2.4mw

## Power requirements


- Two AAA alkaline Batteries
- Power consumption: Less than 40mA
- Battery Life: Two AAA 1.5V, 1200mAh alkaline batteries could be continuously operated as long as 24 hours.

## Environmental requirements

- Operation Temperature: 5°C - 40°C
- Storage Temperature: -25°C - 70°C
- Ambient Humidity: 15% - 93% no condensation in operation; ≤ 93% no condensation in storage/transport
- Atmospheric pressure: 70kPa - 106kPa

## Equipment response time

- As shown in the figure
- Response time of slower average is 8s.

## Classification

- According to the type of protection against electric shock: INTERNALLY POWERED EQUIPMENT;
- According to the degree of protection against electric shock: TYPE BF APPLIED PART (The application part is rubber inside of the Pulse Oximeter);
- According to the degree of protection against ingress of dust and water: IP22
- According to the mode of operation: CONTINUOUS OPERATION

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