

HOMEDICS®

LEADERS IN HEALTH & WELLBEING

✓ Register your product today at
www.homedics.co.uk/product-registration



COMFORT PRO SHAWL

instruction manual

HCM-WRP455H-GB

3 YEAR guarantee

PRODUCT FEATURES:




Fig. 1



Fig. 2



Fig. 3

INSTRUCTIONS FOR USE

1. Insert the jack into the socket on the left side of the wrap (Fig. 1)
2. Plug the adaptor into a AC 220-240v mains outlet.
3. Place massager on your neck and shoulder in the desired position (Fig. 2).
4. To activate the massage action, press the vibration button  on the control pad. The massage will activate on the low setting in steady vibration program (Fig. 3)
5. To change massage intensity, press the vibration button  again
NOTE: there are 3 different intensities, low, medium and high. To flick through them, keep pressing the vibration button.
6. To activate the soothing heat function, press the heat button  on the control pad. The red LED will turn on. To turn off the soothing heat function, simply press the heat button again.
NOTE: You can use the heat without having the vibration.
7. Keep pressing the vibration button until the unit turns off. This will switch off the massage and heat function.
NOTE: There is a 30 minute auto shut-off on this appliance for your safety. This feature should not be considered as a substitute for 'OFF'. Always remember to switch the appliance off when not in use.

MAINTENANCE:

Cleaning

The shawl is wipe clean only using a slightly damp, not wet, cloth. DO NOT immerse in water. DO NOT use petroleum based or flammable cleaning agents because of fire risk. DO NOT use abrasive cleaning agents.

Storage

Switch off the product when not in use. Place it in a safe, dry and dust-free place.

WARNING LABEL ICONS:



Do not wash



NEVER insert pins, metallic fasteners or objects into the appliance or any opening.



Not to be used by children under 3 years of age



Read instructions before use

FKA Brands Ltd

Somerhill Business Park, Tonbridge, Kent TN11 0GP, UK
Customer Support: +44(0) 1732 378557 | support@homedics.co.uk



IB-HCMWRP455H-1018-04