HOMEDICS®

LEADERS IN BEAUTY & WELLBEING

FACIAL BEAUTY ROLLER

WITH MICROCURRENT

3 YEAR GUARANTEE

FAC-BR200-EU2

Register your product today at
www.homedics.co.uk/product-registration

PRODUCT FEATURES:

Fig. 1

1. Massage rollers
2. Solar panel
3. Handle
4. Cleaning cloth
5. Storage bag

HOW IT WORKS

The 360° spinning rollers with microcurrent technology easily slide on the skin, replicating professional kneading massage relieving fatigue.

The 70° V-shape design perfectly fits the facial and body contours for efficient massage. The small, portable, and beautiful design is perfect to add to your daily facial routine and to take with you on the go.

INSTRUCTIONS FOR USE

This device uses a solar panel to generate a microcurrent output to the rollers. There is no need to charge the unit. Either natural or man-made light can be used to power the unit.

NOTE: Make sure you do not cover the solar panel while you are using the device.

You can apply a moisturiser, crème or oil on your skin prior to using the facial roller

Fig. 2

Fig. 3

Fig. 4

Fig. 5

Use it on the skin at different angles to increase the massage efficiency (fig.2).

1. Face: Gently move the Rollers from the jawline to your cheekbone (fig. 3). Afterward, slide the rollers slowly from the corner of your mouth to the ear (fig. 3).
2. Neck: Gently move the Rollers from the shoulder to the neck and jawline (fig. 4).
3. Arms: Gently move the Rollers along the arm (fig. 5).

CAUTION: Do not use this device excessively. It is recommended to massage each area 15-20 times maximum per session.

MAINTENANCE & STORAGE

Clean the device before and after each use to ensure the best performance.

1. Do not immerse in water. Clean the surface of the device with a clean, slightly damp cloth.
2. Never use steel wool, abrasive cleaners or corrosive liquids to clean the device, this may result in malfunction or discolouration.
3. Never use lubricants, detergents, polishing agents or air fresheners on this device.
4. After cleaning, please use a towel to wipe the device.

Storage

Place the appliance in its box or in a safe, dry, cool place. Avoid contact with sharp edges or pointed objects which might damage the surface.

READ ALL INSTRUCTIONS BEFORE USE. SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE.

- This appliance can be used by children aged from 14 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.
- Do not place or store appliance where it can fall or be pulled into a bath or sink. Do not place in or drop into water or other liquid.
- DO NOT reach for an appliance that has fallen in water or other liquids. Keep dry – DO NOT operate in wet or moist conditions.
- NEVER insert pins, metallic fasteners or objects into the appliance or any opening.
- Use this appliance for the intended use as described in this booklet. DO NOT use attachments not recommended by HoMedics.
- NEVER operate the appliance if it is not working properly, if it has been dropped or damaged, or dropped into water. Return it to the HoMedics Service Centre for examination and repair.
- DO NOT operate where aerosol (spray) products are being used or where oxygen is being administered.
- DO NOT attempt to repair the appliance. There are no user serviceable parts. All servicing of this appliance must be performed at an authorised HoMedics Service Centre.
- Please ensure that all hair, clothing and jewellery are kept clear of moving parts of the product at all times.
- If you have any concerns regarding your health, consult a doctor before using this appliance.
- Use of this product should be pleasant and comfortable. Should pain or discomfort result, discontinue use and consult your GP.
- Pregnant women, diabetics and individuals with pacemakers should consult a doctor before using this appliance. Not recommended for use by individuals with sensory deficiencies including diabetic

neuropathy.

- DO NOT use on an infant, invalid or on a sleeping or unconscious person. DO NOT use on insensitive skin or on a person with poor blood circulation.
- This appliance should NEVER be used by any individual suffering from any physical ailment that would limit the user's capacity to operate the controls.
- NEVER use directly on swollen or inflamed areas or skin eruptions.
- Do not use on irritated, sunburned or chapped skin, or on any area with dermatological problems.
- Do not apply over warts, moles or varicose veins.
- Do not use for longer than the recommended time.
- Only gentle force should be exerted against the mechanism in order to eliminate risk of injury.
- Failure to follow the above may result in the risk of fire or injury.

3 YEAR GUARANTEE

FKA Brands Ltd guarantees this product from defect in material and workmanship for a period of 3 years from the date of purchase, except as noted below. This FKA Brands Ltd product guarantee does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorised accessory; alteration to the product; or any other conditions whatsoever that are beyond the control of FKA Brands Ltd. This guarantee is effective only if the product is purchased and operated in the UK / EU. A product that requires modification or adaptation to enable it to operate in any country other than the country for which it was designed, manufactured, approved and / or authorised, or repair of products damaged by these modifications is not covered under this guarantee. FKA Brands Ltd shall not be responsible for any type of incidental, consequential or special damages.

To obtain guarantee service on your product, return the product post-paid to your local service centre along with your dated sales receipt (as proof of purchase). Upon receipt, FKA Brands Ltd will repair or replace, as appropriate, your product and return it to you, post-paid. Guarantee is solely through HoMedics Service Centre. Service of this product by anyone other than the HoMedics Service Centre voids the guarantee. This guarantee does not affect your statutory rights.

For your local HoMedics Service Centre, go to
www.homedics.co.uk/servicecentres

WEEE explanation

This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

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Fig. 1

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Distributed in the UK by: FKA Brands Ltd,
EU Importer: FKA Brands Ltd, 29 Earlsfort Terrace,