instruction manual

**3**YEARguarantee

#### WHAT IS CELLULITE?

Cellulite is a condition in which the skin has a dimpled, lumpy appearance due to underlying fat deposits. It is most noticeable on the buttocks and thighs, but also occurs on the upper arms. By using this product you can help to improve circulation, encourage lymphatic drainage, when combined with eating a healthful, balanced diet and exercising you can help your body to reduce the fat content in cells and reduce the appearance of cellulite.

#### PRODUCT FEATURES

#### CHARGING CABLE

#### Charging

- 1. Plug the jack socket end of the cable into the charging port on the product.
- 2. Plug the USB end of the cable into a standard USB charging plug.
- The LEDs on main unit will light up during charging and they will switch off once fully charged.
- 4. A full charge requires up to 4 hours and will provide up to 45 minutes of use.
- 5. Once fully charged, unplug the adaptor from the socket.

#### INSTRUCTIONS FOR USE

Choose the larger or smaller suction cup, depending on the size of the area to be treated. Attach the suction cup by pushing it firmly onto the product.

- To remove the suction cup, press the small cup release button and pull the cup away from the product.
- Press the LIPO button for 3 seconds to turn on the suction. Press again to increase the suction level. There are 3 levels indicated by the 3 blue LEDs.
- increase the suction level. There are 3 levels indicated by the 3 blue LEDs.

  3. Press the EMS button for 3 seconds to activate the electromagnetic

# stimulation (optional). There are 5 levels indicated by the 5 red LEDs N.B: For the first use, we suggest to start with level 1 for both the EMS and the LIPO.

- There are 2 suctions modes, which can be selected by sliding the button up or down. Suction Mode 1 is intermittent, Suction Mode 2 is continuous.
- 5. Once the settings are as desired, firmly position the cup on the desired body part and glide across the skin very slowly.

N.B: Applying a small amount of a simple baby oil to the skin can assist the gliding effect, but avoid using perfumed oils which may deteriorate the plastic suction cup.

- 6. For optimum results, treatment time is 3 to 5 minutes per area every day.
- 7. When both LIPO and EMS are on, you should feel both suction of your skin and light electrical current. If either becomes uncomfortable, please press the large "Release" button for immediate stop. Adjust the levels to a lower setting before continuing.
- 8. Once done, press the EMS button or LIPO button for 3 seconds to turn the device off. The device includes an auto-shutoff after 10 minutes.

9. See the 'cleaning' section below for advice on how to change the filter

#### **Tips for Use**

- The Cellulite Massager is especially effective, when used regularly and in conjunction with a healthy balanced diet and undertaking of regular exercise.
- During usage, keep moving the product on your body. Just hold firmly the
  unit on the body to avoid dropping it. There is no need to push too hard. If
  the suction power feels uncomfortable and too powerful at any point, push
  "RELEASE" button to release cup from body.
- Some temporary reddening or bruising of the treated area is not uncommon. This should disappear within a few days.

#### **CLEANING:**

- 1. Unplug the appliance before cleaning.
- 2. Remove the cup from the main unit.
- Use the tweezers provided to remove the filter, as shown in the image.
   N.B: we suggest changing the filter once a month or when you start feeling that the suction power is less powerful.
- The suction cup should be wiped clean with a damp cloth after use.
- Never allow water or any other liquids to come into contact with the main unit. Do not immerse in any liquid to clean.
- Never use abrasive cleaners, brushes, glass / furniture polish, paint thinner etc to clean.

#### FKA Brands Ltd

Somerhill Business Park, Tonbridge, Kent TN11 0GP, UK Customer Support: +44(0) 1732 378557 support@homedics.co.uk

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