Register your product today at  
[www.homedicsgroup.com/register](http://www.homedicsgroup.com/register)

# EYE REVIVE LUXE

---

INSTRUCTION MANUAL

3 YEAR Guarantee

EYE-200-EU

**READ ALL INSTRUCTIONS BEFORE USE.  
SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE.**

**IMPORTANT SAFEGUARDS:**

**WHEN USING ELECTRICAL PRODUCTS, ESPECIALLY  
WHEN CHILDREN ARE PRESENT, BASIC SAFETY  
PRECAUTIONS SHOULD ALWAYS BE FOLLOWED,  
INCLUDING THE FOLLOWING:**

- If you have any concerns regarding your health, consult a doctor before using this appliance.
- Only use this appliance for the intended use as described in this booklet.
- This appliance can be used by children aged from 14 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.
- This appliance should NEVER be used by any individual suffering from any physical ailment that would limit the user's capacity to operate the controls.
- Keep away from heated surfaces.
- Individuals with pacemakers should consult a doctor before using this appliance.
- DO NOT use on an infant, invalid or on a sleeping or unconscious person.
- DO NOT use on sensitive, very dry or damaged skin or directly on swollen or inflamed areas or skin eruptions.

If you have circulatory or skin problems consult your doctor before using.

- Do not use on irritated, sunburned or chapped skin, or on any area with dermatological problems.
- If you experience discomfort when using this product, discontinue use. If a more extreme reaction occurs, consult your doctor.
- For external use only.
- This is a non-professional appliance, designed for personal use.
- Only use this appliance around the eyes as described in this booklet.
- Care should be taken when using this device with the red light on. Avoid shining directly into the eye.

**Battery directive**

 This symbol indicates that batteries must not be disposed of in the domestic waste as they contain substances which can be damaging to the environment and health. Please dispose of batteries in designated collection points.

**WEEE explanation**

 This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

## PRODUCT FEATURES:

## USING YOUR EYE MASSAGER

1. Remove any eye make-up and apply a small amount of your favourite eye cream.
  2. Remove the protective cap.
  3. Press and hold the power button to turn the eye massager on.
  4. Toggle through the settings by pressing the power button.
- HOT:** Soothing heat is the perfect setting for relaxing tired eyes and boosting hydration.
- LIGHT:** Combining heat and red light treatments are an effective way to encourage collagen production and reduce fine lines and wrinkles.
- SONIC:** Heat, red light and vibration is the ultimate treatment to refresh and relax your eyes, and lighten the appearance of dark circles.
5. Slide the massage head across your skin as shown, applying light pressure to give a gentle massage. When the correct level of pressure is applied the massager should smoothly glide across your skin, vibrating gently.

6. Press and hold the power button to turn the eye massager off. The eye massager will turn off automatically after 5 minutes of continuous use.
7. Wipe the massage head clean and replace the protective cap after each use.

### Charging your eye massager

1. Plug the micro-USB cable into the input port.
2. Plug the USB end of the cable into a USB power adaptor (not included) or the USB port of a computer.
3. The battery indicator LED will flash to show the massager is being charged.
4. When fully charged the battery indicator LED will turn blue.

## MAINTENANCE:

Regularly check the device for signs of damage. Do not use a damaged device.

### To Store

Store in a safe, dry place, out of children's reach.

### To Clean

Clean the unit after each use with a slightly damp cloth and dry thoroughly before using.

DO NOT use abrasive or chemical cleaners. DO NOT immerse.

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