HOMEDICS®

NOPAIN^{3D}™

Instruction Manual

PRP-200-EU

**2 year
guarantee**

Thank you for choosing HoMedics products. Please read these instructions for use carefully and keep them for further use.

WARNINGS

- Never use this unit near the heart, around the mouth, or on diseased skin.
- Do not use this unit in places with high humidity such as the bathroom.
- User may need to use under the supervision of a physical therapist.
- For external use only.
- Long-term stimulation at the same site may cause skin discomfort.
Consult a dermatologist if skin discomfort persists.
- Do not use the No Pain 3D patch during the first 12 weeks of pregnancy.
- Consult a doctor before you use this device if you receive other medical treatments.
- Do not disassemble or remodel this unit.
- Do not use this unit while sleeping, driving or bathing.
- Prevent any metal object, such as belt buckle or necklace coming into contact with the No Pain 3D patch.
- Supervision and medical consultancy recommended if you are under age 16.

MAINTENANCE & DAMAGE

- If the unit is not functioning properly, stop using immediately.
- If the unit is causing discomfort, stop using immediately.
- Pause the No Pain 3D patch before removing/moving to another area.
- Use only for pain problem diagnosed by professionals.
- Consult your local authorities for proper battery disposal.
- Remove battery from battery chamber if you do not expect to use for a long period of time.

WHAT IS THE NO PAIN 3D PATCH?

The No Pain 3D patch is an innovative one piece self-adhesive electric pain reliever. It uses TENS technology to target specific pain. Built by medical grade silicone, this is ultra-thin, flexible and perfectly matches the contours of contacting surfaces. Advanced circuitry conserves energy by providing 20 x 20 minute treatments per battery life. Above all, it is a safe and easy pain management device that brings people not only a brand new experience, but also a life style.

Read before you use: The No Pain 3D patch should never be used to mask or relieve undiagnosed pain. Before treating any symptoms of pain, it is mandatory to diagnose cause of pain by a doctor.

WHAT IS TENS?


TENS -transcutaneous electrical nerve stimulation - stimulates nerves through the skin. It is recognised as a clinically proven, effective, non-medication method to treat pain from certain causes. It manages pain by stopping the pain gate transmitting pain signals to the brain. It also helps to release endorphins to sooth by stimulating nerves. It is free from side-effects when used properly. The method is scientifically underpinned and medically approved.

APPLICATIONS:


- Back pain
- Joint pain (e.g. knee joint, hip joint, shoulder)
- Neuralgia
- Headache
- Menstrual Pain
- Pain after injuries to musculoskeletal system
- Pain with circulatory problems
- Chronic pain through various causes

ASSEMBLE YOUR NO PAIN 3D PATCH


There are only a few steps of assembly required.


1 Flip open battery cover with a pin with care


2 Insert battery (CR2032 only)

3 Close and lock battery cover


4 Peel the protective film off one side of the gel pads and place adhesive side onto the back of the device. Press firmly.

5 Carefully peel off the other protective film so that the adhesive gel pads are left in place.

HOW TO OPERATE YOUR NO PAIN 3D PATCH

1. Press and hold + to power on the No Pain 3D patch. Smart Contact Detection* is activated with a short and continuous beep
2. Press + again to activate stimulation**
3. Use +/- button to change level of stimulation
4. Automatic power off after 20 minutes
5. Replace battery when the No Pain 3D patch does not turn on.

* Smart Contact Detection automatically detects connection. An audible alarm (constant beeping) will be triggered when there is improper contact between the 2 electrodes.

** When power is on, stimulation level is set to 0 for safety purposes. One more press will activate stimulation mode.


BATTERY REPLACEMENT

Please replace battery when:

1. Power turns off during stimulation and will not turn on again
2. The device will not turn on after a long period of storage.

Note: After Use always place protective film back to electrode after use.

PRODUCT SPECIFICATION

Power: CR2032 Lithium 3V battery

Pulse Rate: 70Hz

Pulse Width: 200 μ s

Output Voltage Max. 65 Vpp, based on 500 Ω load $\pm 10\%$

Treatment Time 20 minutes; auto-off

Pulse Strength 0 ~ 15 stages adjustable

Operation Environment 5 ~ 40°C, 15 ~ 93% Relative Humidity

Storage Environment -25 ~ 70°C, $\leq 93\%$ Relative Humidity

Transportation Environment -25 ~ 70°C, $\leq 93\%$ Relative Humidity

Size 124mm x 41mm x 10mm

PROGRAM MODULATION

Cycle	Mode	Pulse Width	Pulse Rate	On-time
1	Constant	200	70Hz	120s
2	Burst	200	70Hz	120s

If you need longer treatment, please repeat operating procedure.

1. There is no serviceable part content in this product.
2. Replace battery when necessary.

Symbols

The No Pain 3D patch pack and instruction manual are marked with the following symbols:

PRP-200-EU Product code

YYYY-MM Date of manufacture year, month

CE mark with Notified Body registration number (SGS UK: 0120)


Type BF Applied Part symbol ie. insulated patient circuit

Manufacturer complies with the Battery Directive on disposal of used batteries


YYYY-MM Use by date - this relates only to the Self-Adhesive Electrode Gel Pads

Menstrual Pain


Place No Pain 3D patch over the tummy area. For optimum result, move No Pain 3D patch to the other spot after finish stimulation on one spot.

Migraine


Place No Pain 3D patch on the back of the neck. For optimum result, start stimulation at low intensity, then increase gradually until a most suitable intensity achieved.

Note: Do not place pads on the side or front of the neck, or both temples.


Low Back Pain

Place No Pain 3D patch either side at the site of pain. For optimum result, move No Pain 3D patch to the other spot after finish stimulation on one spot.


Back Pain

Place No Pain 3D patch either side of the spine covering the site of pain. For optimum result, start stimulation at low intensity, then increase gradually until a most suitable intensity achieved.

Sciatica

Place No Pain 3D patch either side at the site of pain. For optimum result, move No Pain 3D patch to the other spot after finish stimulation on one spot.

Knee Pain

Place No Pain 3D patch over the top of the knee. Avoid placing directly on the knee cap. For optimum result, start stimulation at low intensity, then increase gradually until a most suitable intensity achieved.

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